



**REBOOTS GO X**



# GET STARTED

## 1 CONNECTING THE CONTROL DEVICE

Attach your Go X control unit to the pants by inserting it into the fastener on the pants and giving it a 90° twist. If you have done everything correctly, the connections on the underside of the Go X control unit will now point towards the waistband. Now insert the plug from the right leg marked with a "R" on the plug into the "R" socket and the plug from the left leg marked with a "L" into the "L" socket.



Please pay strict attention to the correct connection of the left and right plug, otherwise your Go X control unit may become unusable.

## 2 SWITCH ON THE GO X CONTROL UNIT

Switch on the Go X control unit by pressing the "On-Off"-button until you hear a beep. Using the Reboots App with your new Reboots Go X, will give you endless possibilities to create the most personalized Reboots massage ever! So if you are using the App, now is the time to connect your device to your phone.

**DOWNLOAD THE REBOOTS APP HERE:**



App Store



Play Store

## 3 PRESS PLAY AND LEAN BACK

You are ready to go. Press "Play" on the device or choose your preferred program in the app and enjoy your first massage with the Reboots Go X. To get more information on how to use the Reboots Go X, please scan the QR code on the Go X packaging.

**NEED HELP?**



hallo@reboots.de



+49 (0) 6661 748 9106



www.reboots.com